

# Scheduler Module

The Scheduler Module allows you to schedule sessions in advance.

## Dashboard

The dashboard displays a calendar and a list of scheduled and recorded sessions that you are allowed to view. The recorded sessions have a link that allows you to open them.

## Schedule a session

To schedule a session:

1. Open the **Scheduler** tab
2. Double-click a time slot or press the **New schedule** button.
3. Enter a **Session** name in the **Name** field
4. Select the **Location** from the list.
5. Optionally, select other users from the **Participants** list to invite them to the session.
6. Specify the Start and End times.
7. Click **OK**.
8. The scheduled session now appears in the Dashboard of the user that created it and the ones that are invited. [quickstart.scheduler.jpg](#)

---

Revision #3

Created 15 January 2024 11:38:01

Updated 10 March 2024 20:28:52