

# Change the annotator


## aim

To assign the role of session annotator to another user, for instance, a colleague or one of your students.

## prerequisites

- Your license includes the Annotation Module.
- You created the session or received session ownership or you are a Technician or Administrator.

## procedure

1. Open the **Sessions** tab to view the list of all the sessions.
2. Right-click the appropriate session and select **Change annotator**.  
The image shows a small, light blue dialog box with a white background. It has a title bar at the top and a main area with text. The text in the main area reads "Select a user or type unknown". There is a small icon of a person in a circle to the left of the text. The dialog box is positioned over a blurred background of a software interface.
3. From the drop-down list select the name of the person who you want to assign the role of annotator to and click **OK**. By default, the person who recorded the session is the annotator.

---

Revision #2

Created 15 January 2024 10:38:50

Updated 8 March 2024 15:37:17