

Change the annotator


aim

To assign the role of session annotator to another user, for instance, a colleague or one of your students.

prerequisites

- Your license includes the Annotation Module.
- You created the session or received session ownership or you are a Technician or Administrator.

procedure

1. Open the **Sessions** tab to view the list of all the sessions.
2. Right-click the appropriate session and select **Change annotator**.
The image shows a small dialog box with a search bar containing the text 'Search by name or type unknown' and a list of names below it.
3. From the drop-down list select the name of the person who you want to assign the role of annotator to and click **OK**. By default, the person who recorded the session is the annotator.

Revision #2

Created 15 January 2024 11:38:50

Updated 8 March 2024 16:37:17