



# Sessions

The **Sessions** tab contains the list with recorded sessions that you are entitled to review. For Technicians and Administrators, the **Session** tab lists all recorded sessions. For Trainers, only the sessions that they created themselves, received session ownership of, were assigned as annotator for or received an invitation for are listed. For Trainees, the list contains the sessions they received ownership of, were assigned as annotator for or received an invitation for.

 Image not found or type unknown

Double-click a session to open it for review.

See [Review sessions](#)

 Image not found or type unknown

---

Revision #3  
Created 15 January 2024 10:38:46  
Updated 11 March 2024 13:26:38