

# What's new in Viso 12?

## ***Software based licensing fully integrated in Viso***

You can activate or deactivate your software license from the Viso application directly.

## ***Highlight videos***

In the Annotation module, you can now create highlight videos which summarise the most interesting point events of your sessions to be used e.g. in presentations or lectures. To these highlight videos, you may add your markers and remarks as subtitles.

## ***The right to manage recording options moved to Technicians***

Those with a Technician account can now enable multiple location recording, including composite and separate sessions.

## ***Allocate groups to departments***

When creating groups, now the department where the group belongs can also be provided.

## ***Data Transparency Policy***

When a new user logs in, a pop-up window will appear, informing the user about Noldus's Data Transparency Policy.

---

Revision #3

Created 4 March 2024 08:41:19

Updated 28 May 2024 16:38:19 by Niek Wilmink